

NEW

“360°”  
BLADE  
HOLDER

ANY  
ANGLE  
OF USE

SINGLE  
HAND  
LOCKING

- virtually any angle may be created
- blade can be set into a position for any anatomy or angle desired by the surgeon
- particularly advantageous for surgeries like “distal wedge surgeries”
- one-hand opening and closing of locking mechanism
- can be disassembled for cleaning and sterilisation
- registered design and brand TRINOVO<sup>®</sup>

Application videos



Scalpel Blade Holder  
Art.no.: 3303

## Testimonial of the Scalpel Blade Holder - Dr. Ugrinovic & Dr. Sokol



Application images of Dr. Ugrinovic & Dr. Sokol

"The scalpel blade holder is perfect for regions that are difficult to access. The first image shows the access for wisdom tooth removal in the upper jaw. The view of the surgical area is excellent because the hand is far away due to the angle of the grip. The holder also provides a good field of view in the lower jaw, e.g. for wisdom tooth removal. The last two pictures are examples of accesses e.g. when creating a mucosa graft from the palate."

Dr. Sokol – Dr. Ugrinovic & Dr. Sokol, Erlangen (Germany)

**Scalpel Blade Holder**

**Art.no.: 3303**

